



縦磨 kisst (Tatemigakisuto)

This is an eye-opener for you!

The best way to remove plaque is to brush the teeth vertically.

A dentist in Ginza turned his attention to the ordinary toothbrush to solve the fundamental problem with teeth, and spent six long years in clinical trials and experiments to develop this new style of toothbrush.

Please experience the world's first scientifically designed toothbrush created using the best technology by Japanese toothbrush craftsmen.

Different from the usual way of brushing the teeth horizontally adopt **this new idea of brushing vertically** to achieve surprisingly effective plaque removal. Designed to **fit the various alignments of the teeth** (including when straightening). **Reach each tooth easily**. The world's ultimate toothbrush is here now!!

"Stella-arch" (Patent pending)
The new design and shape of the handle allows brushing the backside of the teeth and the molars comfortably and effectively!!
All the teeth are reached easily!

Maximum amount of plaque removal!!

"Unpleasant feeling in the mouth"
That is the indicator!
It's the time to take serious care of the teeth and gums. Quickly experience the dramatic reduction in gingivitis, periodontal disease, or bad breath!!

Q: After eating, I often have food particles stuck between my teeth. I'm worried about this!

>> Food particles stick between the teeth due to recessed gums and wider interdental gaps.

"T-kisst (Tatemigakisuto)" is the toothbrush with the highest usability. Brush the teeth at the point where food particles are stuck between the teeth. No need to use a toothpick, feel good only by brushing.

(T-kisst)BL/PK



Q: I have an unpleasant stickiness in my mouth in the morning, or when I don't talk for a long time!

>> In most cases, this is caused by the growth of bacteria.

"T-kisst" ensures the highest cleaning level for the teeth (plaque removal) compared to any other toothbrush. It removes sugar residue very effectively! Sugar is the nutrient source for the bacteria causing stickiness in the month. Enjoy the exhilarating feeling after brushing the teeth!
Have fun brushing your teeth.

Q: I brush my teeth with an electric toothbrush but have developed a periodontal disease and gingivitis!!

>> You may suffer from hyperaesthesia due to damaged gums caused by applying excessive pressure. Avoid brushing any sore spots as this encourages a vicious circle....

You can remove bacteria stuck in periodontal pockets easily with "T-kisst". Enjoy healthy gums!

We have been waiting for this kind of toothbrush!

Q: Why do my teeth feel smooth for such a long time?

>> This condition lasts a long time because stains sticking to the surface of the teeth, or at the back between the teeth, are removed chemically and completely by the special ion bristles. Enjoy smooth teeth when you get up in the morning. That's the evidence that the oral environment is healthy!

Q: I'm worried about dirty toothbrush bristles.

>> The "T-kisst" uses ion antibacterial bristles manufactured by a special process. After washing in water, the bristles remain clean by preventing the growth of bacteria thanks to the antibacterial properties.

Special direct sales price: 200B each.

● You can brush the teeth by only using water, no need to use toothpaste, a great saving! Thanks to the combination of **the ion antibacterial bristles and natural ceramic apatite components**.

Q: Why is it better to brush the teeth vertically?

>> Plaque between the teeth, or food stuck between the teeth, is removed easily by brushing in the direction of natural tooth growth (vertical)! (Normal brushing for back teeth) The compact head is designed to fit the width of the teeth.

◆ Take care of your teeth to avoid losing them. ◆ 93.87% of 160 users report that "T-kisst" is very effective!

*The new video on how to brush the teeth correctly is now available.

The "T-kisst" toothbrush is designed to make it possible to eat the food you like and keep your teeth for life.

*It is designed for another person to easily brush a child's or an elderly person's teeth when in care.

The person in charge of taking orders in Bangkok: Kikuchi 098-828-5545 [MAIL] dentcosme@outlook.com



Dent+Cosme

Made in Japan [Identification number] 512325794
[HP] <http://www.dentcosme.com/>

